

## Product datasheet for **RC211304L3V**

### **HNRPA3 (HNRNPA3) (NM\_194247) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	HNRPA3 (HNRNPA3) (NM_194247) Human Tagged ORF Clone Lentiviral Particle
Symbol:	HNRPA3
Synonyms:	2610510D13Rik; D10S102; FBRNP; HNRPA3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_194247
ORF Size:	1134 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211304).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_194247.2</a>
RefSeq Size:	5810 bp
RefSeq ORF:	1137 bp
Locus ID:	220988
UniProt ID:	<a href="#">P51991</a>
Cytogenetics:	2q31.2
Protein Pathways:	Spliceosome
MW:	39.6 kDa



[View online »](#)

**Gene Summary:**

Plays a role in cytoplasmic trafficking of RNA. Binds to the cis-acting response element, A2RE. May be involved in pre-mRNA splicing.[UniProtKB/Swiss-Prot Function]