

## Product datasheet for **RC211326L2V**

### NEK9 (NM\_033116) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NEK9 (NM_033116) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NEK9
Synonyms:	APUG; LCCS10; NC; NERCC; NERCC1
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_033116
ORF Size:	2937 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211326).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_033116.3</a>
RefSeq Size:	5560 bp
RefSeq ORF:	2940 bp
Locus ID:	91754
UniProt ID:	<a href="#">Q8TD19</a>
Cytogenetics:	14q24.3
Domains:	RCC1, pkinase, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase


[View online »](#)

**MW:** 107.2 kDa

**Gene Summary:** This gene encodes a member of the NimA (never in mitosis A) family of serine/threonine protein kinases. The encoded protein is activated in mitosis and, in turn, activates other family members during mitosis. This protein also mediates cellular processes that are essential for interphase progression. [provided by RefSeq, Jul 2016]