

## Product datasheet for **RC211340L1V**

### TRIM50 (NM\_178125) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | TRIM50 (NM_178125) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | TRIM50   |
| Synonyms:                 | TRIM50A  |
| Mammalian Cell Selection: | None   |
| Vector:                   | pLenti-C-Myc-DDK (PS100064)  |
| Tag:                      | Myc-DDK  |
| ACCN:                     | NM_178125  |
| ORF Size:                 | 1461 bp  |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC211340).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_178125.2</a> , <a href="#">NP_835226.2</a>  |
| RefSeq Size:              | 2050 bp  |
| RefSeq ORF:               | 1464 bp  |
| Locus ID:                 | 135892   |
| UniProt ID:               | <a href="#">Q86XT4</a>   |
| Cytogenetics:             | 7q11.23  |
| Protein Families:         | Druggable Genome   |
| MW:                       | 54.7 kDa   |


[View online »](#)

**Gene Summary:**

E3 ubiquitin-protein ligase that ubiquitinates Beclin-1/BECN1 in a 'Lys-63'-dependent manner enhancing its binding to ULK1. In turn, promotes starvation-induced autophagy activation. Interacts also with p62/SQSTM1 protein and thereby induces the formation and the autophagy clearance of aggresome-associated polyubiquitinated proteins through HDAC6 interaction.[UniProtKB/Swiss-Prot Function]