

## Product datasheet for **RC211345L3V**

### UTS2D (UTS2B) (NM\_198152) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	UTS2D (UTS2B) (NM_198152) Human Tagged ORF Clone Lentiviral Particle
Symbol:	UTS2D
Synonyms:	U2B; URP; UTS2D
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_198152
ORF Size:	357 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211345).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_198152.3</a>
RefSeq Size:	2408 bp
RefSeq ORF:	360 bp
Locus ID:	257313
UniProt ID:	<a href="#">Q765I0</a>
Cytogenetics:	3q28
Protein Families:	Druggable Genome, Secreted Protein, Transmembrane
MW:	13.7 kDa



[View online »](#)

Gene Summary: Potent vasoconstrictor.[UniProtKB/Swiss-Prot Function]