

Product datasheet for RC211354

LRRC58 (NM_001099678) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	LRRC58 (NM_001099678) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	LRRC58
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC211354 representing NM_001099678 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGCATCGCC

ATGGAGGAGGCCGGAGCAGCGGTGGTCACGGCCGGGGAGGCCGAAGTGAAGTGGTCCCGCCTCAGCGTGT
 CCACCGAGACGCTGGAGTCTGAGCTGGAGGCGCGGGGGAGGAGCGGCGCGGGGCGCGGGAGGCGCTGCT
 GCGGCTGCTGCTGCCTCACAACCGTCTGGTGTGCTGCCACGGGCGCTGGGCAGCGGCTTCCCGCACCTC
 CAGCTGCTGGACGTGAGCGGCAACGCGTTGACCGCGCTCGGGCCGGAGCTGCTCGCTTGCAGCGGCTGC
 GCACGCTGCTGGCCAAGAACAACCGGCTCGGCGGGCCAGTGCCTGCCAAGGGCCTGGCCAGTCGCC
 GCTCTGCCGACGCTCCAGGTGCTCAACCTCAGCGGCAACTGTTCCAGGAGGTGCCTGCCTCGCTCTTA
 GAGCTGCGCGCGCTGCAGACCTGAGCCTGGGCGGCAACCAACTGCAGAGCATCCCGGCTGAGATCGAGA
 ACTTGAGAGTTTAGAGTGTATATCTTGAGGAAATTCATTAAGAAATCCACCGAATTAGGAAA
 TCTGCCTTCTGAATTATTTGGTATTATGTGACAACAAATCCAAAGCATACCTCCTCACTTTACAG
 TTACATTCACCTTCGTTCCCTAAGTCTTCACAATAACTTGCTGACATATCTGCCTCGAGAGATCCTCAACC
 TTATTCATTTGGAAGAGTTGAGTTACGAGGAAATCCATTGGTTGTTGTTTTGTTAGAGATTTAACCTA
 TGATCCTCAACTCTCCTGGAATTAGCTGCACGACCATTAAGATTGAAATATTTCTACACTCCCTAT
 GATCTTCTGAAATCTTCTAGATACTTGGGTTAGCCAGCAATTGCCCAAACCAAGTGTGGAGGAG
 TCTACTTTGACTGCTGTGTGAGACAAATTAATTTGTGGACTTCTGTGGGAAGTATCGCCTCCCACTGAT
 GCACTACTTGTGTTACACAGAATGTTCTTCCCCTTGCAAGTTCTGCCTCTCACAGTCCACTTCCAGAGT
 GAATCTGACTCAGAAGATGAAGTAGTGTGCTGCACGCAGAATGCAGAAAGTTCTTCTTGGT

AGCGGACCGACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC211354 representing NM_001099678
 Red=Cloning site Green=Tags(s)


MEEAGAAVVTAGEAELNWSRLSVSTETLESELEARGEERRGAREALLRLLLPHNRLVSLPRALGSGFPHL
 QLLDVSGNALTALGPELLALRGLRTLLAKNNRLGGPSALPKGLAQSPLCRSLQVLNLSGNCFQEVASLL
 ELRALQTL SLGGNQLQSIPAEIENLQSLECLYLGGNF IKEIPPELGNLPSLNYLVLCDNKIQSIPPQLSQ
 LHSLRSLSLHNNLLTYLPREILNLIHLEELSLRGNPLVVRFVRDLTYDPPTLLELAARTIKIRNISYTPY
 DLPGNLLRYLGSASNCNPNKCGGVYFDCCVRQIKFVDFCGKYRLPLMHYLCSPCCSSPCSSASHSSTSQS
 ESDSEDEASVAARRMQKVLGG

SGP^{TRTRRL}EQKLI^{SEEDLAANDILDYKDDDDK}V

Chromatograms: https://cdn.origene.com/chromatograms/mk8025_d08.zip

Restriction Sites: SgfI-RsrII

Cloning Scheme:



ACCN: NM_001099678

ORF Size: 1113 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001099678.2

RefSeq Size: 7683 bp

RefSeq ORF: 1116 bp

Locus ID: 116064


UniProt ID: Q96CX6

Cytogenetics: 3q13.33


Protein Families: Druggable Genome

MW: 40.4 kDa

Product images:



Circular map for RC211354



Western blot validation of overexpression lysate (Cat# [LY420490]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC211354 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).