

Product datasheet for **RC211354L4V**

LRRC58 (NM_001099678) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LRRC58 (NM_001099678) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LRRC58
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001099678
ORF Size:	1113 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211354).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001099678.1
RefSeq Size:	7683 bp
RefSeq ORF:	1116 bp
Locus ID:	116064
UniProt ID:	Q96CX6
Cytogenetics:	3q13.33
Protein Families:	Druggable Genome
MW:	40.4 kDa



[View online »](#)