

## Product datasheet for RC211357L2V

### MYBPHL (NM\_001010985) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	MYBPHL (NM_001010985) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MYBPHL
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001010985
ORF Size:	1062 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211357).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001010985.1</a>
RefSeq Size:	1372 bp
RefSeq ORF:	1065 bp
Locus ID:	343263
UniProt ID:	<a href="#">A2RUH7</a>
Cytogenetics:	1p13.3
MW:	38.7 kDa
Gene Summary:	This gene encodes a protein with two immunoglobulin superfamily domains and a fibronectin 3 domain. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2012]


[View online »](#)