

Product datasheet for **RC211372**

NRG1 (NM_013962) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	NRG1 (NM_013962) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	NRG1
Synonyms:	ARIA; GGF; GGF2; HGL; HRG; HRG1; HRGA; MST131; MSTP131; NDF; NRG1-IT2; SMDF
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>RC211372 representing NM_013962
 Red=Cloning site Blue=ORF Green=Tags(s)

CTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCCGGCGC
 GCC

ATGAGATGGCGACGCGCCCCGCGCGCTCCGGGCGTCCCGCCCCGGGCCAGCGCCCCGGCTCCGCCG
 CCGCTCGTCGCCGCGCTGCCGCTGCTGCCACTACTGCTGCTGCTGGGACCGCGGCCCTGGCGCCGG
 GCGCGCGGCCGCAACGAGGCGGCTCCCGCGGGGCTCGGTGTGCTACTCGTCCCCGCCAGCGTGGGA
 TCGGTGCAGGAGCTAGCTCAGCGCGCCGCGGTGGTGTGATCGAGGAAAGGTGCACCCGAGCGGGCAGC
 AGGGGGCACTCGACAGGAAGGCGCGCGCGGGCGAGGCAGGGGCGTGGGGCGCGATCGCGAGCC
 GCCAGCCGCGGGCCACGGGCGCTGGGGCCGCCCGGAGGAGCCGCTGCTCGCCGCAACGGGACCGTG
 CCCTCTTGGCCACCGCCCCGTGCCAGCGCCGCGAGCCGGGAGGAGGCCCTATCTGGTGAAGG
 TGCACCAGGTGTGGGCGGTGAAAGCCGGGGCTTGAAGAAGGACTCGCTGCTCACCGTGCCTGGGGAC
 CTGGGGCCACCCGCCTTCCCCTCTGCGGGAGGCTCAAGGAGGACAGCAGGTACATCTTCTTCATGGAG
 CCCGACGCCAACAGCACCAGCCGCGCGCGGCCCTTCCGAGCCTTTTCCCCCTCTGGAGACGGGCC
 GGAACCTCAAGAAGGAGGTGAGCCGGTGTGTGCAAGCGGTGCGCCTTGCTCCCAATTGAAAGAGAT
 GAAAAGCCAGGAATCGGCTGCAGGTTCCAACTAGTCCTTCGGTGTGAAACCAGTTCTGAATACTCCTCT
 CTCAGATTCAGTGGTTCAGAAATGGGAATGAATTGAATCGAAAAACAAACCACAAAATATCAAGATAC
 AAAAAAGCCAGGGAAGTCAGAACTTCGCATTAACAAAGCATCACTGGCTGATTCTGGAGAGTATATGTG
 CAAAGTGATCAGCAAATAGGAAATGACAGTGCCTCTGCCAATATCACCATCGTGAATCAAACGCTACA
 TCTACATCCACTGGGACAAGCCATCTTGTAAAATGTGCGGAGAAGGAGAAAACCTTCTGTGTGAATG
 GAGGGGAGTGCTTCATGGTGAAGACCTTCAAACCCCTCGAGATACTTGTGCAAGTGCCCAAATGAGTT
 TACTGGTATCGCTGCCAAAACACGTAATGCCAGCTTCTACAGTACGTCCACTCCCTTCTGTCTCTG
 CCTGAA

ACGCGTACGCGCGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RC211372 representing NM_013962
 Red=Cloning site Green=Tags(s)

MRWRRAPRRSGRPGPRAQRPGSAARSSPPLPLPLLLLLLGTAAALPGAAAGNEAAPAGASVCYSSPPSVG
 SVQELAQRAAVVIEGKVHPQRRQGGALDRKAAAAAGEAGAWGGDREPPAAGPRALGPPAEPEPLAANGTV
 PSWPTAPVPSAGEPGEEAPYLKVVHVVAVKAGGLKKDSLLTVRLGTWGHAPAFPSGRLKEDSRVYFFME
 PDANSTRAPAAFRAFPPLETGRNLKKEVSRVLCRRCALPPQLKEMKSQESAAGSKLVLCRSETSSSEYSS
 LRFKWFKNGNELNRKNKPQNIKIQQKPGKSELRINKASLADSGEYMCKVISKLGNDASANITIVESNAT
 STSTTGTSHLVKCAEKEKTFVNGGECFMVKDLNPSRYLCKCPNEFTGDRQNYMASFYSTSTPFLSL
 PE


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms:

https://cdn.origene.com/chromatograms/mg4516_h08.zip

Restriction Sites:

Ascl-MluI

Cloning Scheme:


ACCN: NM_013962

ORF Size: 1266 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_013962.1](#), [NP_039256.1](#)

RefSeq Size: 1986 bp

RefSeq ORF: 1269 bp

Locus ID: 3084

UniProt ID: [Q02297](#)

Cytogenetics: 8p12


Protein Families: Druggable Genome, Secreted Protein, Transcription Factors, Transmembrane

Protein Pathways: ErbB signaling pathway

MW: 45 kDa

Gene Summary: The protein encoded by this gene is a membrane glycoprotein that mediates cell-cell signaling and plays a critical role in the growth and development of multiple organ systems. An extraordinary variety of different isoforms are produced from this gene through alternative promoter usage and splicing. These isoforms are expressed in a tissue-specific manner and differ significantly in their structure, and are classified as types I, II, III, IV, V and VI. Dysregulation of this gene has been linked to diseases such as cancer, schizophrenia, and bipolar disorder (BPD). [provided by RefSeq, Apr 2016]

Product images:



Circular map for RC211372