

Product datasheet for RC211395

MDFIC (NM_199072) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	MDFIC (NM_199072) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	MDFIC
Synonyms:	HIC; MDFIC1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC211395 representing NM_199072. Blue=ORF Red=Cloning site Green=Tag(s)

```
GCTCGTTTGTAGTGAACCGTCAGAATTTGTAAATACGACTACTATAGGGCGGCCGGGAATTCGTGCGACTG
GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
GTGCGCGGAGTCAGAGCCGCCACCGCTGCCGAGTTGCCGCCACTGCGGCGTCTGGGCTGAGCCGGAGG
GAGGCGGGAGGACGCGCAGGGCGGCCGCCGCGCTCGTCAGGCCACCGGGCGAAAATGCGGCCGCTGC
CGGAGGCTCGCTAACTTTCCGGGGCGGAAGAGGAGGAGGAGGAGGAGGAAGGGCTTGGAGCGACTACG
GGGGATGCGGAGAAGCAGTCAGTTCCTGCACCCAGCACCTCACAGCCCTTCTCCGTGCGCCCTGCC
GGGCGGCGAGCTAGGCGGCAGCGCGCGCGGGCTCGGCGGAGCGGCCATGTCCGGCGGGCGAA
GCCCTCGCTCCCGGGCCCGTGGGGCCGAGCGCGTGGCCGAGGCGGGCGGCCAGCTGGGCTCCACA
GCCCAGGGAAAATGTGATAAAGACAATACTGAGAAAGATATAACTCAAGTACCAATAGCCACTTCACA
CATGGAGAGATGCAAGACCAGTCCATTTGGGAAATCCTTCGGATGGTGAATCATTAGAACCCAACT
CAGCGCTTGCCTCAGCTTCAGACTTCAGCCAGGTGCCAAGTGGTGGAGAAATAGGCAAGATAAAGAAC
GGCCACACAGGTCTGAGCAATGGAATGGAATTCACCACGGGGCCAAACCGGATCCGCAGATAATCGC
AACTTTCAGCACCTGTTTCTCAAAAAATGCATAGAAAAATTCAGTCCAGCTTGTCTGTAACAGCGAT
ATCAGTAAGAAGAGCAAAGTAAATGCTGTCTTTCCAAAAGACAGGCTTTCACCTGAAGATTGTTGT
GTCCACTGTATCCTGGCTTGTCTGTTCTGCGAATTCCTGACCTTTCACACATTGCTCTGGGACAAGCG
TCATGTGGCATCTGCACCTCAGAAGCCTGCTGCTGTTGCTGTGGTGACGAGATGGGGATGATTGTAA
TGCCCTTGTGATATGACTGTGGCATCATGGATGCCTGTTGTGAATCATCAGACTGCTTGGAAATCTGT
ATGGAATGCTGTGGAATTTGTTTCCTTCA
ACGCGTACGCGGCCGCTCGAGCAGAAAATCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGAT
TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC
```



[View online »](#)


Protein Sequence: >Peptide sequence encoded by RC211395
 Blue=ORF Red=Cloning site Green=Tag(s)

VRGVRAATAAAVAATAASGLSRREAGGRAGAAAAVVRPPGRKCGRCRRLANFPGRKRRRRRRKGLGATT
 GGCGEAVSSLHPAPHSPSSVRPAGRRARRQRRGAGSAERPMSGAGEALAPGPVGPQRVAEAGGGQLGST
 AQGKCDKDNTKIDITQATNSHFTHGEMQDQSIWGNPSDGELIRTQPQLPQLQTSAQVPSGEEIGKIKN
 GHTGLSNGNGIHHGAKHGSADNRKLSAPVVSQKMHRKIQSSLVNSDI SKKSKVNAVFSQKTGSSPEDCC
 VHCILACLFCFLTL CNIVL GQASCGICTSEACCCCGDEMGGDCNCPDMDCGIMDACCESSDCLCIC
 MECCGICFPS
 TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8025_d05.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_199072

ORF Size: 1065 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_199072.5](#)

RefSeq Size: 4602 bp

RefSeq ORF: 1068 bp

Locus ID: 29969


UniProt ID: [Q9P1T7](#)

Cytogenetics: 7q31.1-q31.2

MW: 37 kDa

Gene Summary: This gene product is a member of a family of proteins characterized by a specific cysteine-rich C-terminal domain, which is involved in transcriptional regulation of viral genome expression. Alternative translation initiation from an upstream non-AUG (GUG), and an in-frame, downstream AUG codon, results in the production of two isoforms, p40 and p32, respectively, which have different subcellular localization; p32 is mainly found in the cytoplasm, whereas p40 is targeted to the nucleolus. Both isoforms have transcriptional regulatory activity that is attributable to the cysteine-rich C-terminal domain. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]

Product images:



Circular map for RC211395