

Product datasheet for RC211408L4

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com

OriGene Technologies, Inc.

nttps://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

FOXA2 (NM_153675) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: FOXA2 (NM 153675) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: FOXA2

Synonyms: HNF-3-beta; HNF3B; TCF3B

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC211408).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_153675

ORF Size: 1371 bp





FOXA2 (NM_153675) Human Tagged Lenti ORF Clone - RC211408L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 153675.1</u>

 RefSeq Size:
 2230 bp

 RefSeq ORF:
 1374 bp

 Locus ID:
 3170

 UniProt ID:
 Q9Y261

Cytogenetics: 20p11.21

Protein Families: Embryonic stem cells, ES Cell Differentiation/IPS, Induced pluripotent stem cells, Transcription

Factors

Protein Pathways: Maturity onset diabetes of the young

MW: 48.1 kDa


Gene Summary: This gene encodes a member of the forkhead class of DNA-binding proteins. These

hepatocyte nuclear factors are transcriptional activators for liver-specific genes such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. This gene has been linked to sporadic cases of maturity-onset diabetes of the young. Transcript variants encoding different isoforms have been identified for this gene.


[provided by RefSeq, Oct 2008]



Product images:



Circular map for RC211408L4



Double digestion of RC211408L4 using Sgfl and Mlul