

Product datasheet for **RC211468L2V**

DCAF17 (NM_025000) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | DCAF17 (NM_025000) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | DCAF17 |
| Synonyms: | C2orf37; C20orf37 |
| Mammalian Cell Selection: | None |
| Vector: | pLenti-C-mGFP (PS100071) |
| Tag: | mGFP |
| ACCN: | NM_025000 |
| ORF Size: | 1560 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC211468). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_025000.2 |
| RefSeq Size: | 5744 bp |
| RefSeq ORF: | 1563 bp |
| Locus ID: | 80067 |
| UniProt ID: | Q5H9S7 |
| Cytogenetics: | 2q31.1 |
| MW: | 58.6 kDa |



[View online »](#)

Gene Summary:

This gene encodes a nuclear transmembrane protein that associates with cullin 4A/damaged DNA binding protein 1 ubiquitin ligase complex. Mutations in this gene are associated with Woodhouse-Sakati syndrome. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Sep 2009]