

Product datasheet for **RC211475A1V**

Human SIAH1 (NM_001006610) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human SIAH1 (NM_001006610) AAV Particle
Tag:	Myc-DDK
Symbol:	SIAH1
Synonyms:	BURHAS; SIAH1A
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC211475 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACGGGAAAGGCTACTCCACCTTCTCTGTACTCCTGGAGGGGAGTCTTGTTACATGTTTACCAGCGG
CCAGGACAAGGAAGAGAAAAGAAATGAGCCGTCAGACTGCTACAGCATTACCTACCGGTACCTCGAAGTG
TCCACCATCCCAGAGGGTGCCCTGCCCTGACTGGCACAACCTGCATCCAACAATGACTTGCGGAGTCTTTTT
GAGTGTCCAGTCTGCTTTGACTATGTGTTACCGCCCATTTCTCAATGTCAGAGTGCCCATCTTGTTTGA
GCAACTGTCGCCCAAAGCTCACATGTTGTCCAACCTGCCGGGGCCCTTTGGGATCCATTGCAACTTGGC
TATGGAGAAAGTGGCTAATTTCAGTACTTTTCCCTGTAAATATGCGTCTTCTGGATGTGAAATAACTCTG
CCACACACAGAAAAGCAGACCATGAAGAGCTCTGTGAGTTTAGGCCTTATTCCTGTCCGTGCCCTGGTG
CTTCCTGTAAATGGCAAGGCTCTCTGGATGCTGTAATGCCCATCTGATGCATCAGCATAAGTCCATTAC
AACCTACAGGGAGAGGATATAGTTTTCTTGCTACAGACATTAATCTTCCTGGTGCTGTTGACTGGGTG
ATGATGCAGTCTGTTTTGGCTTTCACCTCATGTTAGTCTTAGAGAAACAGGAAAAATACGATGGTCACC
AGCAGTTCTTCGCAATCGTACAGCTGATAGGAACACGCAAGCAAGCTGAAAAATTTGCTTACCGACTTGA
GCTAAATGGTCATAGGCGACGATTGACTTGGGAAGCGACTCCTCGATCTATTTCATGAAGGAATTGCAACA
GCCATTATGAATAGCGACTGTCTAGTCTTTGACACCAGCATTGCACAGCTTTTGCAGAAAATGGCAATT
TAGGCATCAATGTAACATTTTCATGTGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC211475 protein sequence Red=Cloning site Green=Tags(s) MTGKATPPSLYSWRGVLTCLPAARTRKRKEMSRQTATALPTGTSKCPPSQRPALTGTTASNNDLASLF ECPVCFDYVLPPIQCQSGHLVCSNCRPKLTCCPTCRGPLGSIRNLAMEKVANSVLPCKYASSGCEITL PHTEKADHEELCEFRPYSCPCPGASCKWQGS�DAVMPHLMHQHSITTLQGEDIVFLATDINLPGAVDWV MMQSCFGFHFMLVLEKQEKYDGHQQFFAIVQLIGTRKQAEINFAYRLELNGHRRRLTWEATPRSIHEGIAT AIMNSDCLVFDTIAQLFAENGLGINVTISM TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001006610
ORF Size:	939 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001006610.1</u>
RefSeq Size:	2393 bp
RefSeq ORF:	942 bp
Locus ID:	6477
UniProt ID:	<u>Q8IUQ4</u>
Cytogenetics:	16q12.1
MW:	34.6 kDa