

Product datasheet for **RC211484L3V**

NEK4 (NM_003157) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	NEK4 (NM_003157) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NEK4
Synonyms:	NRK2; pp12301; STK2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003157
ORF Size:	2523 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211484).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003157.1
RefSeq Size:	3698 bp
RefSeq ORF:	2526 bp
Locus ID:	6787
UniProt ID:	A5YM70
Cytogenetics:	3p21.1
Domains:	ppkinase, TyrKc, S_TKc
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

MW: 94.4 kDa

Gene Summary: The protein encoded by this gene is a serine/threonine protein kinase required for normal entry into replicative senescence. The encoded protein also is involved in cell cycle arrest in response to double-stranded DNA damage. Finally, this protein plays a role in maintaining cilium integrity, and defects in this gene have been associated with ciliopathies. [provided by RefSeq, Jan 2017]