

Product datasheet for RC211524

BRCC45 (BRE) (NM_199192) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BRCC45 (BRE) (NM_199192) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	BRCC45
Synonyms:	BRCC4; BRCC45; BRE
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC211524 representing NM_199192 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCCAGAAGTGGCCTTGAACCGAATATCTCCAATGCTCTCCCCTTTCATATCTAGCGTGGTCCGGA
ATGAAAAGTGGGACTGGATGCTACAACTGTTTGAGGATAACTGACTTAAAATCTGGCTGCACATCATT
GACTCTGGGCCCAACTGTGACCGATTTAACTGCACATACCATATGCTGGAGAGACATTAAGTGGGAT
ATCATTTCATGCCCAATACCCAGAAGTGCCTCCGATTTTATCTTTGGAGAAGATGCTGAATTCCTGC
CAGACCCCTCAGCTTTCGAGAATCTTGCTCCTGGAATCCTTCAAATCCTGAATGTCTTACTTGTGGT
GAAGGAATTGTGCAACAATATACCAATTCCAATGTAGCCGCTCCGGGAGAGCTCCCGCTCATGTTT
GAATACCAGACATTACTGGAGGACACAGTATGGAGAGAATGGAAATTTATGCTGGGAAAAAACA
ACTGGACTGGTGAATTTTCAGCTCGTTTCTTTGAAAGCTGCCCGTAGATTTTCAGCAATATCCCCACATA
CCTTCTCAAGGATGTAATGAAGACCCTGGAGAAGATGTGGCCCTCCTCTCTGTTAGTTTGGAGCACT
GAAGCCACCCAGGTGTACCCAGCTGTACTTGTACCTCGAATTGAGCATGCACCTGGAGGCTCCTCAG
CTCTTCATATCCAGCTTTTCCAGGAGGAGGATGTCTCATTGATTACGTTCTCAAGTATGCCACCTGCT
CACCAACAAGGTGCAGTACGTGATTCAAGGGTATCACAAAAGAAGAGAGTATATTGCTGCTTTTCTCAGT
CACTTTGGCACAGGTGCTGGAATATGATGCAGAAGGCTTACAAAACACTCACTGCTGCTGATGTGGA
AAGATTTTGTCTTGTACACATTGACCTGCCTCTGTTTTTCCCTCGAGACCAGCCAACCTCTCACATT
TCAGTCCGTTTATCACTTTACCAACAGTGGACAGCTTTACTCCAGGCCAAAAAATTATCCGTACAGC
CCCAGATGGGATGAAATGAAATGGCCAAAAGAGCAAAAAGAGAGCAACAGAGATGGGAGGAATCCA
GCTCTGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC211524 representing NM_199192
Red=Cloning site Green=Tags(s)

MSPEVALNRISPMLSPFISSVVRNGKVGLDATNCLRITDLKSGCTSLTPGPNCDRFKLHIPPYAGETLKWD
 IIFNAQYPELPPDFIFGEDAEFLPDPALQNLASWNPSNPECLLLVVKELVQYHQFQCSRLRESSRLMF
 EYQTLLEEPQYGENMEIYAGKKNWTGEFSARFLKLPVDFSNIPTYLLKDVNEDPGEDVALLSVSFEDT
 EATQVYPKLYLSPRIEHALGGSSALHIPAFPPGGCLIDYVPQVCHLLTNKVQYVIQGYHKRREYIAAFLS
 HFGTGVVEYDAEGFTKLTLMLMKWDFCFLVHIDLPLFFPRDQPTLTFQSVYHFTNSGQLYSQAQKNYPYS
 PRWDGNEMAKRAKRESNRDGEESSA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk8055_f10.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



* The last codon before the Stop codon of the ORF

ACCN: NM_199192

ORF Size: 1128 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_199192.2](#), [NP_954662.1](#)

RefSeq Size: 1791 bp

RefSeq ORF: 1131 bp

Locus ID: 9577

UniProt ID: [Q9NXR7](#)


Cytogenetics: 2p23.2

Protein Families: Druggable Genome

MW: 42.5 kDa

Gene Summary: This gene encodes an anti-apoptotic, death receptor-associated protein that interacts with tumor necrosis factor-receptor-1. The encoded protein acts as an adapter in several protein complexes, including the BRCA1-A complex and the BRISC complex. The BRCA1-A complex possesses ubiquitinase activity and targets sites of double strand DNA breaks, while the BRISC complex exhibits deubiquitinase activity and is involved in mitotic spindle assembly. This gene is upregulated in several types of cancer. [provided by RefSeq, Jun 2016]

Product images:



Circular map for RC211524