

Product datasheet for **RC211528L1V**

GABA A Receptor gamma 2 (GABRG2) (NM_000816) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GABA A Receptor gamma 2 (GABRG2) (NM_000816) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GABA A Receptor gamma 2
Synonyms:	CAE2; DEE74; ECA2; EIEE74; FEB8; GEFSP3
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_000816
ORF Size:	1401 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211528).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000816.2
RefSeq Size:	2370 bp
RefSeq ORF:	1404 bp
Locus ID:	2566
UniProt ID:	P18507
Cytogenetics:	5q34
Domains:	Neur_chan_memb, Neur_chan_LBD



[View online »](#)

Protein Families:	Druggable Genome, Ion Channels: Cys-loop Receptors, Transmembrane
Protein Pathways:	Neuroactive ligand-receptor interaction
MW:	54.16 kDa
Gene Summary:	This gene encodes a gamma-aminobutyric acid (GABA) receptor. GABA is the major inhibitory neurotransmitter in the mammalian brain, where it acts at GABA-A receptors, which are ligand-gated chloride channels. GABA-A receptors are pentameric, consisting of proteins from several subunit classes: alpha, beta, gamma, delta and rho. Mutations in this gene have been associated with epilepsy and febrile seizures. Multiple transcript variants encoding different isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]