

Product datasheet for RC211547L3

ROCK1 (NM_005406) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: ROCK1 (NM_005406) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: ROCK1

Synonyms: P160ROCK; ROCK-I

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC211547).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





 $[\]ensuremath{^*}$ The last codon before the Stop codon of the ORF.




OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn



Plasmid Map:



ACCN: NM_005406 **ORF Size:** 4062 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts

of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at customercom or by

calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Cytogenetics:

ROCK1 (NM_005406) Human Tagged Lenti ORF Clone - RC211547L3

Reconstitution Method:

- 1. Centrifuge at 5,000xg for 5min.
- 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
- 3. Close the tube and incubate for 10 minutes at room temperature.
- 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: <u>NM 005406.2</u>

 RefSeq Size:
 6650 bp

 RefSeq ORF:
 4065 bp

 Locus ID:
 6093

 UniProt ID:
 Q13464

Domains: pkinase, HR1, PH

Protein Families: Druggable Genome, Protein Kinase

18q11.1

Protein Pathways: Axon guidance, Chemokine signaling pathway, Focal adhesion, Leukocyte transendothelial

migration, Pathogenic Escherichia coli infection, Regulation of actin cytoskeleton, TGF-beta

signaling pathway, Vascular smooth muscle contraction, Wnt signaling pathway

MW: 158.6 kDa

Gene Summary: This gene encodes a protein serine/threonine kinase that is activated when bound to the

GTP-bound form of Rho. The small GTPase Rho regulates formation of focal adhesions and stress fibers of fibroblasts, as well as adhesion and aggregation of platelets and lymphocytes by shuttling between the inactive GDP-bound form and the active GTP-bound form. Rho is also essential in cytokinesis and plays a role in transcriptional activation by serum response factor. This protein, a downstream effector of Rho, phosphorylates and activates LIM kinase,

which in turn, phosphorylates cofilin, inhibiting its actin-depolymerizing activity. A

pseudogene, related to this gene, is also located on chromosome 18. [provided by RefSeq,

Aug 2015]