

Product datasheet for **RC211566L4V**

TRAF3IP1 (NM_015650) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TRAF3IP1 (NM_015650) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TRAF3IP1
Synonyms:	FAP116; IFT54; MIP-T3; MIPT3; SLSN9
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_015650
ORF Size:	2073 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211566).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_015650.2
RefSeq Size:	2568 bp
RefSeq ORF:	2076 bp
Locus ID:	26146
UniProt ID:	Q8TDR0
Cytogenetics:	2q37.3
Protein Families:	Druggable Genome
MW:	78.5 kDa


[View online »](#)

Gene Summary:

The protein encoded by this gene interacts with TNF receptor-associated factor 3, tethering it to cytoskeletal microtubules. The encoded protein is also an inhibitor of the innate type I IFN response. Defects in this gene are a cause of Senior-Loken syndrome 9. [provided by RefSeq, Mar 2017]