

Product datasheet for **RC211610L4V**

ARHGAP17 (NM_018054) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ARHGAP17 (NM_018054) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARHGAP17
Synonyms:	MST066; MST110; MSTP038; MSTP066; MSTP110; NADRIN; PP367; PP4534; RICH-1; RICH1; WBP15
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018054
ORF Size:	2409 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211610).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018054.4
RefSeq Size:	3246 bp
RefSeq ORF:	2412 bp
Locus ID:	55114
UniProt ID:	Q68EM7
Cytogenetics:	16p12.1
Domains:	RhoGAP, BAR



[View online »](#)

Protein Families: Druggable Genome

MW: 87.4 kDa

Gene Summary: RICH1 is a GTPase-activating protein (GAP). GAPs stimulate the intrinsic GTP hydrolysis of small G proteins, such as RHOA (MIM 165390), RAC1 (MIM 602048), and CDC42 (MIM 116952). [supplied by OMIM, Apr 2004]