

Product datasheet for **RC211625L3V**

TOX2 (NM_001098798) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TOX2 (NM_001098798) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TOX2
Synonyms:	C20orf100; dj495O3.1; dj1108D11.2; GCX-1; GCX1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001098798
ORF Size:	1464 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211625).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001098798.1
RefSeq Size:	2428 bp
RefSeq ORF:	1467 bp
Locus ID:	84969
UniProt ID:	Q96NM4
Cytogenetics:	20q13.12
Protein Families:	Transcription Factors
MW:	51.4 kDa



[View online »](#)

Gene Summary: Putative transcriptional activator involved in the hypothalamo-pituitary-gonadal system.
[UniProtKB/Swiss-Prot Function]