

## Product datasheet for **RC211635L3V**

### BCR (NM\_004327) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	BCR (NM_004327) Human Tagged ORF Clone Lentiviral Particle
Symbol:	BCR
Synonyms:	ALL; BCR1; CML; D22S11; D22S662; PHL
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004327
ORF Size:	3813 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211635).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_004327.3</a>
RefSeq Size:	4739 bp
RefSeq ORF:	3816 bp
Locus ID:	613
UniProt ID:	<a href="#">P11274</a>
Cytogenetics:	22q11.23
Domains:	C2, RhoGAP, RhoGEF, PH
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

<b>Protein Pathways:</b>	Chronic myeloid leukemia, Pathways in cancer
<b>MW:</b>	142.6 kDa
<b>Gene Summary:</b>	<p>A reciprocal translocation between chromosomes 22 and 9 produces the Philadelphia chromosome, which is often found in patients with chronic myelogenous leukemia. The chromosome 22 breakpoint for this translocation is located within the BCR gene. The translocation produces a fusion protein which is encoded by sequence from both BCR and ABL, the gene at the chromosome 9 breakpoint. Although the BCR-ABL fusion protein has been extensively studied, the function of the normal BCR gene product is not clear. The unregulated tyrosine kinase activity of BCR-ABL1 contributes to the immortality of leukaemic cells. The BCR protein has serine/threonine kinase activity and is a GTPase-activating protein for p21rac and other kinases. Two transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Jan 2020]</p>