

Product datasheet for **RC211643L3V**

CSHL1 (NM_022580) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CSHL1 (NM_022580) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CSHL1
Synonyms:	CS-5; CSHP1; CSL; GHB4; hCS-L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022580
ORF Size:	417 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211643).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022580.1
RefSeq Size:	588 bp
RefSeq ORF:	420 bp
Locus ID:	1444
UniProt ID:	Q14406
Cytogenetics:	17q23.3
Domains:	hormone
Protein Families:	Secreted Protein



[View online »](#)

MW: 15.8 kDa

Gene Summary: The protein encoded by this gene is a member of the somatotropin/prolactin family of hormones which play an important role in growth control. The gene, along with four other related genes, is located at the growth hormone locus on chromosome 17 where they are interspersed in the same transcriptional orientation; an arrangement which is thought to have evolved by a series of gene duplications. Although the five genes share a remarkably high degree of sequence identity, they are expressed selectively in different tissues. This particular family member is expressed in placental villi, although it was originally thought to be a pseudogene. In fact, alternative splicing suggests that the majority of the transcripts would be unable to express a secreted protein. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]