

## Product datasheet for **RC211649L4V**

### C1QTNF8 (NM\_207419) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	C1QTNF8 (NM_207419) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C1QTNF8
Synonyms:	CTRP8; UNQ5829
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_207419
ORF Size:	756 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211649).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_207419.2</a>
RefSeq Size:	2974 bp
RefSeq ORF:	759 bp
Locus ID:	390664
UniProt ID:	<a href="#">P60827</a>
Cytogenetics:	16p13.3
Protein Families:	Secreted Protein
MW:	27.5 kDa



[View online »](#)

**Gene Summary:** May play a role as ligand of RXFP1.[UniProtKB/Swiss-Prot Function]