

## Product datasheet for **RC211651L4V**

### Histidase (HAL) (NM\_002108) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Histidase (HAL) (NM_002108) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Histidase
Synonyms:	HIS; HSTD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_002108
ORF Size:	1971 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211651).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_002108.2</a>
RefSeq Size:	3927 bp
RefSeq ORF:	1974 bp
Locus ID:	3034
UniProt ID:	<a href="#">P42357</a>
Cytogenetics:	12q23.1
Protein Families:	Druggable Genome
Protein Pathways:	Histidine metabolism, Metabolic pathways, Nitrogen metabolism



[View online »](#)

**MW:** 72.7 kDa

**Gene Summary:** Histidine ammonia-lyase is a cytosolic enzyme catalyzing the first reaction in histidine catabolism, the nonoxidative deamination of L-histidine to trans-urocanic acid. Histidine ammonia-lyase defects cause histidinemia which is characterized by increased histidine and histamine and decreased urocanic acid in body fluids. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Apr 2012]