

Product datasheet for RC211656A1V

Human OCLM (NM_022375) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human OCLM (NM_022375) AAV Particle
Tag: Myc-DDK
Symbol: OCLM
Synonyms: TISR
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC211656 representing NM_022375
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGGGATGTATCCACCATTATTGTTAAAGATTACCTTAGCAGACACATTAGTATACTCTTCTATTTAA
 AAATCCTTTATAAGAGTGGTATTATATGGCTATCCTGGTATTCTTTCATATTGTTAGTACTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC211656 representing NM_022375
 Red=Cloning site Green=Tags(s)

MGMYPPLLLKIYLSRHISILFYLKILYKSGIWLWYSFILLVL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_022375
ORF Size: 132 bp
Buffer: PBS with 0.001% Pluronic F68


[View online »](#)

Stability: AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.

RefSeq: [NM_022375.2](#), [NP_071770.1](#)

RefSeq Size: 884 bp

RefSeq ORF: 135 bp

Locus ID: 10896

Cytogenetics: 1q31.1

MW: 5.1 kDa