

Product datasheet for **RC211665L3V**

KSR2 (NM_173598) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | KSR2 (NM_173598) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | KSR2 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-Myc-DDK-P2A-Puro (PS100092) |
| Tag: | Myc-DDK |
| ACCN: | NM_173598 |
| ORF Size: | 2763 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC211665). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_173598.4 , NP_775869.3 |
| RefSeq Size: | 17008 bp |
| RefSeq ORF: | 2853 bp |
| Locus ID: | 283455 |
| UniProt ID: | Q6VAB6 |
| Cytogenetics: | 12q24.22-q24.23 |
| Protein Families: | Druggable Genome, Protein Kinase |
| MW: | 104.1 kDa |



[View online »](#)

Gene Summary:

Location-regulated scaffold connecting MEK to RAF. Has very low protein kinase activity and can phosphorylate MAP2K1 at several Ser and Thr residues with very low efficiency (in vitro). Interaction with BRAF enhances KSR2-mediated phosphorylation of MAP2K1 (in vitro). Blocks MAP3K8 kinase activity and MAP3K8-mediated signaling. Acts as a negative regulator of MAP3K3-mediated activation of ERK, JNK and NF-kappa-B pathways, inhibiting MAP3K3-mediated interleukin-8 production.[UniProtKB/Swiss-Prot Function]