

Product datasheet for **RC211699L3V**

TMEM219 (NM_001083613) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM219 (NM_001083613) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM219
Synonyms:	IGFBP-3R
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001083613
ORF Size:	720 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211699).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001083613.1
RefSeq Size:	1019 bp
RefSeq ORF:	723 bp
Locus ID:	124446
UniProt ID:	Q86XT9
Cytogenetics:	16p11.2
Protein Families:	Druggable Genome, Transmembrane
MW:	25.5 kDa



[View online »](#)

Gene Summary: Cell death receptor specific for IGFBP3, may mediate caspase-8-dependent apoptosis upon ligand binding.[UniProtKB/Swiss-Prot Function]