

Product datasheet for **RC211711L3V**

LYNX1 (SLURP2) (NM_177458) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LYNX1 (SLURP2) (NM_177458) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLURP2
Synonyms:	SLURP-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_177458
ORF Size:	291 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211711).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_177458.1
RefSeq Size:	589 bp
RefSeq ORF:	294 bp
Locus ID:	432355
UniProt ID:	P0DP57
Cytogenetics:	8q24.3
MW:	10 kDa



[View online »](#)

Gene Summary:

This gene encodes a novel, secreted member of the Ly6/uPAR (LU) superfamily of proteins containing the unique three-finger LU domain. This gene is mainly expressed in epithelial cells, including skin and keratinocytes, and is up-regulated in psoriatic skin lesions, suggesting its involvement in the pathophysiology of psoriasis. Alternatively spliced transcript variants have been found for this gene. Read-through transcription from the neighboring upstream gene (LYNX1) generates naturally-occurring transcripts (LYNX1-SLURP2) that encode a fusion protein comprised of sequence sharing identity with each individual gene product. [provided by RefSeq, Sep 2017]