

## Product datasheet for RC211711L4

### LYNX1 (SLURP2) (NM\_177458) Human Tagged Lenti ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	LYNX1 (SLURP2) (NM_177458) Human Tagged Lenti ORF Clone
Tag:	mGFP
Symbol:	SLURP2
Synonyms:	SLURP-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
E. coli Selection:	Chloramphenicol (34 ug/mL)
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211711).
Restriction Sites:	SgfI-MluI
Cloning Scheme:	

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF.


ACCN:	NM_177458
ORF Size:	291 bp



[View online »](#)

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_177458.1</a>
<b>RefSeq Size:</b>	589 bp
<b>RefSeq ORF:</b>	294 bp
<b>Locus ID:</b>	432355
<b>UniProt ID:</b>	<a href="#">P0DP57</a>
<b>Cytogenetics:</b>	8q24.3
<b>MW:</b>	10 kDa
<b>Gene Summary:</b>	This gene encodes a novel, secreted member of the Ly6/uPAR (LU) superfamily of proteins containing the unique three-finger LU domain. This gene is mainly expressed in epithelial cells, including skin and keratinocytes, and is up-regulated in psoriatic skin lesions, suggesting its involvement in the pathophysiology of psoriasis. Alternatively spliced transcript variants have been found for this gene. Read-through transcription from the neighboring upstream gene (LYNX1) generates naturally-occurring transcripts (LYNX1-SLURP2) that encode a fusion protein comprised of sequence sharing identity with each individual gene product. [provided by RefSeq, Sep 2017]

Product images:



Circular map for RC211711L4