

Product datasheet for RC211715

ATP5F1D (NM_001687) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: ATP5F1D (NM_001687) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: ATP5F1D
Synonyms: ATP5D; MC5DN5
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC211715 representing NM_001687
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGCTGCCCGCCGCGCTGCTCCGCCGCCGGGACTTGCCCGCCTCGTCCGCCACGCCCGTGCCTATGCCG
 AGGCCGCCCGCCCGGCTGCCGCTCTGCCCAACCAGATGTCTTACCTTCGCTCTCCACGCA
 GGTGTTCTTCAACGGTGCCAACGTCCGGCAGGTGGACGTGCCACGCTGACCGGAGCCTTCGGCATCCTG
 GCGGCCACGTGCCACGCTGCAGGTCTGCGGCCGGGCTGGTGGTGCATGCAGAGGACGGCACCA
 CCTCAAATACTTTGTGAGCAGCGTTCATCGCAGTGAACGCCGACTTTCGGTGCAGTTGTTGCCGA
 AGAGGCCGTGACGCTGGACATGTTGGACCTGGGGCAGCCAAGGCAAACCTGGAGAAGGCCAGGCGGAG
 CTGGTGGGACAGCTGACGAGGCCACGGGCGAGATCCAGATCCGAATCGAGGCCAACGAGGCCCTGG
 TGAAGGCCCTGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC211715 representing NM_001687
 Red=Cloning site Green=Tags(s)

MLPAALLRRPGLGRLVRHARAYAEAAAAPAAASGPNQMSFTFASPTQVFFNGANVRQVDVPTLTGAFGIL
 AAHVPTLQVLRPGLVVVHAEDGTTSKYFVSSGSIAVNADSSVQLLAAEAVTLDMLDLGAANKLEKAQAE
 LVGTADEATRAEIQIRIEANEALVKALE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Chromatograms: https://cdn.origene.com/chromatograms/mk6764_b06.zip



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_001687

ORF Size: 504 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_001687.5](#)


RefSeq Size: 1005 bp

RefSeq ORF: 507 bp


Locus ID: 513
UniProt ID: [P30049](#)
Cytogenetics: 19p13.3
Domains: ATP-synt_DE
Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease
MW: 17.49 kDa
Gene Summary:

This gene encodes a subunit of mitochondrial ATP synthase. Mitochondrial ATP synthase catalyzes ATP synthesis, utilizing an electrochemical gradient of protons across the inner membrane during oxidative phosphorylation. ATP synthase is composed of two linked multi-subunit complexes: the soluble catalytic core, F1, and the membrane-spanning component, Fo, comprising the proton channel. The catalytic portion of mitochondrial ATP synthase consists of 5 different subunits (alpha, beta, gamma, delta, and epsilon) assembled with a stoichiometry of 3 alpha, 3 beta, and a single representative of the other 3. The proton channel consists of three main subunits (a, b, c). This gene encodes the delta subunit of the catalytic core. Alternatively spliced transcript variants encoding the same isoform have been identified. [provided by RefSeq, Jul 2008]


Product images:



Circular map for RC211715



Western blot validation of overexpression lysate (Cat# [LY419790]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC211715 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified ATP5F1D protein (Cat# [TP311715]). The protein was produced from HEK293T cells transfected with ATP5F1D cDNA clone (Cat# RC211715) using MegaTran 2.0 (Cat# [TT210002]).