

## Product datasheet for RC211735L4

#### **OriGene Technologies, Inc.** 9620 Medical Center Drive, Ste 200

Rockville, MD 20850, US
Phone: +1-888-267-4436
https://www.origene.com
techsupport@origene.com
EU: info-de@origene.com
CN: techsupport@origene.cn

### XRCC3 (NM\_001100118) Human Tagged Lenti ORF Clone

#### **Product data:**

**Product Type:** Expression Plasmids

Product Name: XRCC3 (NM\_001100118) Human Tagged Lenti ORF Clone

Tag:mGFPSymbol:XRCC3Synonyms:CMM6

Mammalian Cell Puromycin

Selection:

**Vector:** pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

**ORF Nucleotide** The ORF insert of this clone is exactly the same as(RC211735).

Sequence:

Restriction Sites: Sgfl-Mlul

**Cloning Scheme:** 





<sup>\*</sup> The last codon before the Stop codon of the ORF

ACCN: NM\_001100118

ORF Size: 1038 bp





#### XRCC3 (NM\_001100118) Human Tagged Lenti ORF Clone - RC211735L4

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001100118.1</u>, <u>NP 001093588.1</u>

RefSeq Size: 2563 bp
RefSeq ORF: 1041 bp
Locus ID: 7517
UniProt ID: 043542

**Cytogenetics:** 14q32.33

**Protein Families:** Druggable Genome

**Protein Pathways:** Homologous recombination

**MW:** 37.9 kDa

**Gene Summary:** This gene encodes a member of the RecA/Rad51-related protein family that participates in

homologous recombination to maintain chromosome stability and repair DNA damage. This gene functionally complements Chinese hamster irs1SF, a repair-deficient mutant that


exhibits hypersensitivity to a number of different DNA-damaging agents and is

chromosomally unstable. A rare microsatellite polymorphism in this gene is associated with cancer in patients of varying radiosensitivity. Alternatively spliced transcript variants encoding

the same protein have been identified. [provided by RefSeq, Jul 2008]



# **Product images:**



Circular map for RC211735L4