

Product datasheet for RC211759

CSHL1 (NM_022581) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: CSHL1 (NM_022581) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: CSHL1
Synonyms: CS-5; CSHP1; CSL; GHB4; hCS-L
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC211759 representing NM_022581
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTGCAGGCTCCCGGACGTCCTGCTCCTGGCTTTTGGCCTGCTCTGCCTGCCCTGGCTTCAAGAGG
 CTGGTGCCGTCCAAACCGTTCCCTTATCCAGGCTTTTTAAAGAGGCTATGCTCCAAGCCATCGCGCACA
 CCAGCTGGCCATTGACACCTACCAGGAGTTTATAAGCTCTTGGGAATGGACTCTATTCCGACATCCTCC
 AACATGGAGGAAACGCAGCAGAAATCCAACCTAGAGCTGCTCCACATCTCCCTGCTGCTCATCGAGTCG
 GGCTGGAGCCCGTCCGGTTCCTCAGGAGTACCTCACCAACAACCTGGTGTATGACACCTCGGACAGCGA
 TGAATACCTCCTAAAGGACCTAGAGGAAGGCATCAAATGCTGATGGGGAGGCTGGAAGACGGCAGC
 CACCTGACTGGGACAGCCCTCAAGCAGACCTACAGCAAGTTTGACACAACTCGCACAACCATGACGCAC
 TGCTCAAGAACTACGGGCTGCTCCACTGCTTCAGGAAGGACATGGACAAGGTCGAGACATTCTCGCAT
 GGTGCAGTGCCGCTCTGTGGAGGCGAGCTGTGGCTTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC211759 representing NM_022581
 Red=Cloning site Green=Tags(s)

MAAGSRTSLLAFALLCLPWLQEAGAVQTVPLSRLFKEAMLQAHRAHQLAIDTYQEFISSWGMSIPTSS
 NMEETQQKSNLELLHISLLIESRLEPVRFLRSTFTNNLVYDTSDDYHLLKDLEEGIQMLMGRLEDGS
 HLTGQTLKQTYSKFDTNSHNHDALLKNYGLLHCFRKMMDKVETFLRMVQCRSVEGSCGF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_022581

ORF Size: 597 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

- Reconstitution Method:
1. Centrifuge at 5,000xg for 5min.
 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
 3. Close the tube and incubate for 10 minutes at room temperature.
 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_022581.3](#)

RefSeq Size: 768 bp

RefSeq ORF: 600 bp

Locus ID: 1444

UniProt ID: [Q14406](#)


Cytogenetics: 17q23.3

Protein Families: Secreted Protein

MW: 19.9 kDa

Gene Summary: The protein encoded by this gene is a member of the somatotropin/prolactin family of hormones which play an important role in growth control. The gene, along with four other related genes, is located at the growth hormone locus on chromosome 17 where they are interspersed in the same transcriptional orientation; an arrangement which is thought to have evolved by a series of gene duplications. Although the five genes share a remarkably high degree of sequence identity, they are expressed selectively in different tissues. This particular family member is expressed in placental villi, although it was originally thought to be a pseudogene. In fact, alternative splicing suggests that the majority of the transcripts would be unable to express a secreted protein. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]

Product images:



Circular map for RC211759