

Product datasheet for **RC211796L2V**

VLDL Receptor (VLDLR) (NM_003383) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	VLDL Receptor (VLDLR) (NM_003383) Human Tagged ORF Clone Lentiviral Particle
Symbol:	VLDL Receptor
Synonyms:	CAMRQ1; CARMQ1; CHRMQ1; VLDL-R; VLDLRCH
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_003383
ORF Size:	2619 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211796).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_003383.3
RefSeq Size:	3646 bp
RefSeq ORF:	2622 bp
Locus ID:	7436
UniProt ID:	P98155
Cytogenetics:	9p24.2
Domains:	Idl_recept_b, EGF_CA, Idl_recept_a, EGF
Protein Families:	Druggable Genome, Transmembrane



[View online »](#)

MW: 96.1 kDa

Gene Summary: The low density lipoprotein receptor (LDLR) gene family consists of cell surface proteins involved in receptor-mediated endocytosis of specific ligands. This gene encodes a lipoprotein receptor that is a member of the LDLR family and plays important roles in VLDL-triglyceride metabolism and the reelin signaling pathway. Mutations in this gene cause VLDLR-associated cerebellar hypoplasia. Alternative splicing generates multiple transcript variants encoding distinct isoforms for this gene. [provided by RefSeq, Aug 2009]