

Product datasheet for RC211797L3

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Oxytocin Receptor (OXTR) (NM_000916) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Oxytocin Receptor (OXTR) (NM_000916) Human Tagged Lenti ORF Clone

Tag: Myc-DDK

Symbol: Oxytocin Receptor

Synonyms: OT-R; OTR

Mammalian Cell Puromycin

Selection:

Vector:

pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone


Sequence:

The ORF insert of this clone is exactly the same as(RC211797).

Restriction Sites: Sgfl-Rsrll

Cloning Scheme:





st The last codon before the Stop codon of the ORF.

ACCN: NM_000916

ORF Size: 1176 bp



Oxytocin Receptor (OXTR) (NM_000916) Human Tagged Lenti ORF Clone - RC211797L3

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000916.3</u>, <u>NP 000907.2</u>

 RefSeq Size:
 4361 bp

 RefSeq ORF:
 1170 bp

 Locus ID:
 5021

 UniProt ID:
 P30559

Domains: 7tm 1

Cytogenetics:

Protein Families: Druggable Genome, GPCR, Transmembrane

3p25.3

Protein Pathways: Calcium signaling pathway, Neuroactive ligand-receptor interaction

MW: 42.6 kDa

Gene Summary: The protein encoded by this gene belongs to the G-protein coupled receptor family and acts

as a receptor for oxytocin. Its activity is mediated by G proteins which activate a

phosphatidylinositol-calcium second messenger system. The oxytocin-oxytocin receptor system plays an important role in the uterus during parturition. [provided by RefSeq, Jul

2008]