

Product datasheet for RC211824L4

GRK7 (NM_139209) Human Tagged Lenti ORF Clone

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

Product data:

Product Type: Expression Plasmids

Product Name: GRK7 (NM_139209) Human Tagged Lenti ORF Clone

Tag: mGFP
Symbol: GRK7
Synonyms: GPRK7

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC211824).

Sequence:

Restriction Sites: Sgfl-Rsrll

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_139209

ORF Size: 1659 bp





GRK7 (NM_139209) Human Tagged Lenti ORF Clone - RC211824L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 139209.2</u>, <u>NP 631948.1</u>

 RefSeq Size:
 1746 bp

 RefSeq ORF:
 1662 bp

 Locus ID:
 131890

 UniProt ID:
 Q8WTQ7

Cytogenetics: 3q23

Protein Families: Druggable Genome, Protein Kinase

Protein Pathways: Chemokine signaling pathway, Endocytosis

MW: 62 kDa


Gene Summary: This gene encodes a member of the guanine nucleotide-binding protein (G protein)-coupled

receptor kinase subfamily of the Ser/Thr protein kinase family. It is specifically expressed in the retina and the encoded protein has been shown to phosphorylate cone opsins and

initiate their deactivation. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC211824L4