

## Product datasheet for **RC211837L4V**

### SCO2 (NM\_005138) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

|                           |  |
|---------------------------|--|
| Product Type:             | Lentiviral Particles   |
| Product Name:             | SCO2 (NM_005138) Human Tagged ORF Clone Lentiviral Particle  |
| Symbol:                   | SCO2   |
| Synonyms:                 | CEMCOX1; ECGF1; Gliostatin; MC4DN2; MYP6; PD-ECGF; SCO1L; TdRPase; TP; TYMP  |
| Mammalian Cell Selection: | Puromycin  |
| Vector:                   | pLenti-C-mGFP-P2A-Puro (PS100093)  |
| Tag:                      | mGFP   |
| ACCN:                     | NM_005138  |
| ORF Size:                 | 798 bp   |
| ORF Nucleotide Sequence:  | The ORF insert of this clone is exactly the same as(RC211837).   |
| OTI Disclaimer:           | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a> |
| OTI Annotation:           | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.   |
| RefSeq:                   | <a href="#">NM_005138.1</a>  |
| RefSeq Size:              | 1005 bp  |
| RefSeq ORF:               | 801 bp   |
| Locus ID:                 | 9997   |
| UniProt ID:               | <a href="#">O43819</a>   |
| Cytogenetics:             | 22q13.33   |
| Domains:                  | SCO1-SenC  |
| Protein Families:         | Druggable Genome   |



[View online »](#)

**MW:** 29.8 kDa

**Gene Summary:** Cytochrome c oxidase (COX) catalyzes the transfer of electrons from cytochrome c to molecular oxygen, which helps to maintain the proton gradient across the inner mitochondrial membrane that is necessary for aerobic ATP production. Human COX is a multimeric protein complex that requires several assembly factors; this gene encodes one of the COX assembly factors. The encoded protein is a metallochaperone that is involved in the biogenesis of cytochrome c oxidase subunit II. Mutations in this gene are associated with fatal infantile encephalocardiomyopathy and myopia 6. [provided by RefSeq, Oct 2014]