

Product datasheet for RC211875L1

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Clusterin (CLU) (NM_001831) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Clusterin (CLU) (NM_001831) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: Clusterin

Synonyms: AAG4; APO-J; APOJ; CLI; CLU1; CLU2; KUB1; NA1/NA2; SGP-2; SGP2; SP-40; TRPM-2; TRPM2

Mammalian Cell None

Selection:


Vector:pLenti-C-Myc-DDK (PS100064)E. coli Selection:Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC211875).

Sequence:

Restriction Sites: Sgfl-Rsrll

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_001831

ORF Size: 1503 bp





Clusterin (CLU) (NM_001831) Human Tagged Lenti ORF Clone - RC211875L1

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001831.2, NP 001822.2</u>

 RefSeq Size:
 2859 bp

 RefSeq ORF:
 1350 bp

 Locus ID:
 1191

 UniProt ID:
 P10909

Cytogenetics: 8p21.1

Domains:

Protein Families: Druggable Genome, Secreted Protein

CLa

MW: 57.7 kDa


Gene Summary: The protein encoded by this gene is a secreted chaperone that can under some stress

conditions also be found in the cell cytosol. It has been suggested to be involved in several basic biological events such as cell death, tumor progression, and neurodegenerative disorders. Alternate splicing results in both coding and non-coding variants. [provided by


RefSeq, May 2011]



Product images:



Circular map for RC211875L1



Double digestion of RC211875L1 using Sgfl and Rsrll