

## Product datasheet for **RC211884L4V**

### **MKS1 (NM\_017777) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	MKS1 (NM_017777) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MKS1
Synonyms:	BBS13; JBTS28; MES; MKS; POC12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017777
ORF Size:	1677 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211884).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017777.2</a>
RefSeq Size:	2402 bp
RefSeq ORF:	1680 bp
Locus ID:	54903
UniProt ID:	<a href="#">Q9NXB0</a>
Cytogenetics:	17q22
MW:	64.3 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene localizes to the basal body and is required for formation of the primary cilium in ciliated epithelial cells. Mutations in this gene result in Meckel syndrome type 1 and in Bardet-Biedl syndrome type 13. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2009]