

## Product datasheet for RC211911L1V

### PTCD1 (NM\_015545) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PTCD1 (NM_015545) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PTCD1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_015545
ORF Size:	2100 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211911).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_015545.2</a>
RefSeq Size:	3079 bp
RefSeq ORF:	2103 bp
Locus ID:	26024
UniProt ID:	<a href="#">O75127</a>
Cytogenetics:	7q22.1
Domains:	PPR
MW:	78.7 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a mitochondrial protein that binds leucine tRNAs and other mitochondrial RNAs and plays a role in the regulation of translation. Increased expression of this gene results in decreased mitochondrial leucine tRNA levels. Naturally occurring read-through transcription exists between upstream ATP5J2 (ATP synthase, H<sup>+</sup> transporting, mitochondrial Fo complex, subunit F2) and this gene. [provided by RefSeq, Aug 2015]