

Product datasheet for RC211924L4

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

OriGene Technologies, Inc.

MEIS2 (NM_172315) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: MEIS2 (NM_172315) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: MEIS2

Synonyms: CPCMR; HsT18361; MRG1

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC211924).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_172315

ORF Size: 1164 bp





MEIS2 (NM_172315) Human Tagged Lenti ORF Clone - RC211924L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 172315.1</u>

RefSeq Size: 2775 bp
RefSeq ORF: 1167 bp
Locus ID: 4212

UniProt ID: O14770

Cytogenetics: 15q14

Protein Families: Transcription Factors

MW: 42.1 kDa


Gene Summary: This gene encodes a homeobox protein belonging to the TALE ('three amino acid loop

extension') family of homeodomain-containing proteins. TALE homeobox proteins are highly conserved transcription regulators, and several members have been shown to be essential contributors to developmental programs. Multiple transcript variants encoding distinct

isoforms have been described for this gene. [provided by RefSeq, Jul 2008]



Product images:



Circular map for RC211924L4