

Product datasheet for **RC211974A1V**

Human CAPSL (NM_144647) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human CAPSL (NM_144647) AAV Particle
Tag: Myc-DDK
Symbol: CAPSL
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC211974 ORF sequence
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCGATCCAGGCCAAGAAAAAGCTCACACGGCCACCAACCCATTGAAAGACTCCGACTGCAGTGCC
TGGCCAGGGGCTCTGCTGGGATCAAAGGACTTGGCAGAGTGTTTAGAATTATGGATGACGATAATAATCG
AACCCCTGATTTTAAAGAATTTATGAAAGGGTTAAATGATTATGCTGTGGTCATGGAAAAAGAAGAGGTG
GAAGAAGCTTTCCGGAGGTTTGATAAAGATGGAATGGAACAATAGACTTCAATGAATTTCTTCTCACAT
TAAGACCTCCAATGTCCAGAGCCAGAAAAGAGTAATCATGCAAGCTTTTAGAAAGTTAGACAAGACTGG
AGATGGTGTTATAACAATCGAAGACCTTCGTGAAGTATATAATGCAAAACACCACCCAAAGTACCAGAAT
GGGGAATGGAGTGAGGAACAAGTATTTAGGAAATTTCTGGATAACTTTGATTACCCTATGACAAAGATG
GATTGGTGACCCCTGAGGAGTTCATGAAGTACTATGCAGGTGTGAGCGCATCCATTGACACTGATGTGTA
CTTCATCATCATGATGAGAACCCTGGAAGCTT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC211974 protein sequence
Red=Cloning site Green=Tags(s)

MAIQAKKKLTATNPIERLRLQCLARGSAGIKGLGRVFRIMDDNNRTLDFKEFMKGLNDYAVVMEKEEV
EELFRFDKDGNGTIDFNEFLLTRPPMSRARKEVIMQAFRLDKTGDGVITIEDLREVYNKHHPKYQN
GEWSEEQVFRKFLDNFSPYDKDGLVTPEEFMNYAGVSASIDTDVYFIIMMRTAWKL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2



[View online »](#)

ACCN:	NM_144647
ORF Size:	594 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_144647.3, NP_653248.2</u>
RefSeq Size:	1019 bp
RefSeq ORF:	627 bp
Locus ID:	133690
UniProt ID:	<u>Q8WWF8</u>
Cytogenetics:	5p13.2
MW:	23.1 kDa