

## Product datasheet for **RC211996L3V**

### **PAX2 (NM\_003990) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	PAX2 (NM_003990) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PAX2
Synonyms:	FSGS7; PAPRS
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_003990
ORF Size:	1293 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC211996).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_003990.2</a>
RefSeq Size:	4257 bp
RefSeq ORF:	1299 bp
Locus ID:	5076
UniProt ID:	<a href="#">Q02962</a>
Cytogenetics:	10q24.31
Protein Families:	Druggable Genome
MW:	46.1 kDa



[View online »](#)

**Gene Summary:**

PAX2 encodes paired box gene 2, one of many human homologues of the *Drosophila melanogaster* gene *prd*. The central feature of this transcription factor gene family is the conserved DNA-binding paired box domain. PAX2 is believed to be a target of transcriptional suppression by the tumor suppressor gene WT1. Mutations within PAX2 have been shown to result in optic nerve colobomas and renal hypoplasia. Alternative splicing of this gene results in multiple transcript variants. [provided by RefSeq, Dec 2014]