

## Product datasheet for **RC212019L2V**

### WIPF1 (NM\_001077269) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	WIPF1 (NM_001077269) Human Tagged ORF Clone Lentiviral Particle
Symbol:	WIPF1
Synonyms:	PRPL-2; WAS2; WASPIP; WIP
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001077269
ORF Size:	1509 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212019).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001077269.1</a> , <a href="#">NP_001070737.1</a>
RefSeq Size:	4664 bp
RefSeq ORF:	1512 bp
Locus ID:	7456
UniProt ID:	<a href="#">O43516</a>
Cytogenetics:	2q31.1
MW:	51.1 kDa


[View online »](#)

**Gene Summary:**

This gene encodes a protein that plays an important role in the organization of the actin cytoskeleton. The encoded protein binds to a region of Wiskott-Aldrich syndrome protein that is frequently mutated in Wiskott-Aldrich syndrome, an X-linked recessive disorder. Impairment of the interaction between these two proteins may contribute to the disease. Two transcript variants encoding the same protein have been identified for this gene. [provided by RefSeq, Jul 2008]