

Product datasheet for **RC212040L3V**

EPB41 (NM_004437) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EPB41 (NM_004437) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EPB41
Synonyms:	4.1R; EL1; HE
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004437
ORF Size:	1764 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212040).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004437.3
RefSeq Size:	5823 bp
RefSeq ORF:	1767 bp
Locus ID:	2035
UniProt ID:	P11171
Cytogenetics:	1p35.3
Domains:	B41, SAB
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Tight junction

MW: 66.4 kDa

Gene Summary: The protein encoded by this gene, together with spectrin and actin, constitute the red cell membrane cytoskeletal network. This complex plays a critical role in erythrocyte shape and deformability. Mutations in this gene are associated with type 1 elliptocytosis (EL1). Alternatively spliced transcript variants encoding different isoforms have been described for this gene.[provided by RefSeq, Oct 2009]