

Product datasheet for **RC212044A1V**

Human Eotaxin 3 (CCL26) (NM_006072) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Eotaxin 3 (CCL26) (NM_006072) AAV Particle
Tag: Myc-DDK
Symbol: Eotaxin 3
Synonyms: IMAC; MIP-4a; MIP-4alpha; SCYA26; TSC-1
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC212044 representing NM_006072
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGATGGGCTCTCCTTGGCCTCTGCTGTGCTCCTGGCCTCCCTCCTGAGTCTCCACCTTGGAACTGCCA
CACGTGGGAGTGACATATCCAAGACCTGCTGCTTCCAATACAGCCACAAGCCCTTCCCTGGACCTGGGT
GCGAAGCTATGAATTCACCAAGTAAACAGCTGCTCCAGCGGGCTGTGATATTCCTACCAAAAGAGGCAAG
AAAGTCTGTACCCATCCAAGGAAAAAATGGGTGCAAAAATACATTCTTTACTGAAAACCTCCGAAACAAT
TG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC212044 representing NM_006072
Red=Cloning site Green=Tags(s)

MMGLSLASAVLLASLLSLHLGTATRGSDISKTCFQYSHKPLPWTWVRSYEFTSNSCSQRAVIFTTKRGK
KVCTHPRKKWVQKYISLLKTPKQL

TRTRPLE**QKL**ISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_006072
ORF Size: 282 bp



[View online »](#)

Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006072.4, NP_006063.1</u>
RefSeq Size:	562 bp
RefSeq ORF:	285 bp
Locus ID:	10344
UniProt ID:	<u>Q9Y258</u>
Cytogenetics:	7q11.23
MW:	10.65 kDa