

Product datasheet for **RC212069A1V**

Human HRK (NM_003806) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human HRK (NM_003806) AAV Particle
Tag: Myc-DDK
Symbol: HRK
Synonyms: DP5; HARAKIRI
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC212069 representing NM_003806
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGTGCCCGTGCCCCCTGCACCGCGGCCGCGCCCCCGGCCGTGTGCGCCTGCAGCGCGGGTCGCCTGG
GGCTGCGCTCGTCCGCCGCGCAGCTACCGCGCGCCGGCTCAAGGCGCTAGGCGACGAGCTGCACGCG
CACCATGTGGCGGCGCGCGCGGAGCCGAGGGCGCGCGCCGCGCGCTCCACCTACTGGCCT
TGGCTGTGCGCGGCCGCGCAGGTGGCGCGCTGGCGGCCGTGGCTGCTCGGCAGGCGGAATTG

ACGCGTACGCGGCGGCTCGAGCAGAACTCATCTCTGAAGAGGATCTGGCAGCAATGATATCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC212069 representing NM_003806
Red=Cloning site Green=Tags(s)

MCPCPLHRGRGPPAVCACSAGRLGLRSSAAQLTAARLKALGDELHQRTMWRRRARSRRAPAGALPTYWP
WLCAAAQVAALAAWLLGRRNL*

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Human
Serotype: AAV-2
ACCN: NM_003806
ORF Size: 273 bp
Buffer: PBS with 0.001% Pluronic F68



[View online »](#)

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_003806.1</u>
RefSeq Size:	716 bp
RefSeq ORF:	276 bp
Locus ID:	8739
UniProt ID:	<u>O00198</u>
Cytogenetics:	12q24.22
MW:	9.7 kDa