

Product datasheet for RC212112L4

PAX2 (NM_000278) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: PAX2 (NM_000278) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: PAX2

Synonyms: FSGS7; PAPRS

Mammalian Cell Puromycin

Selection:

Vector: pLenti-C-mGFP-P2A-Puro (PS100093)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC212112).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF

ACCN: NM_000278

ORF Size: 1179 bp



OriGene Technologies, Inc. 9620 Medical Center Drive, Ste 200

CN: techsupport@origene.cn

Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com



PAX2 (NM_000278) Human Tagged Lenti ORF Clone - RC212112L4

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 000278.2</u>

RefSeq Size: 4207 bp
RefSeq ORF: 1185 bp
Locus ID: 5076
UniProt ID: Q02962

 UniProt ID:
 Q02962

 Cytogenetics:
 10q24.31

Protein Families: Druggable Genome

MW: 41.9 kDa


Gene Summary: PAX2 encodes paired box gene 2, one of many human homologues of the Drosophila

melanogaster gene prd. The central feature of this transcription factor gene family is the conserved DNA-binding paired box domain. PAX2 is believed to be a target of transcriptional supression by the tumor suppressor gene WT1. Mutations within PAX2 have been shown to result in optic nerve colobomas and renal hypoplasia. Alternative splicing of this gene results

in multiple transcript variants. [provided by RefSeq, Dec 2014]



Product images:



Circular map for RC212112L4