

Product datasheet for RC212135

Claudin18 (CLDN18) (NM_016369) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Claudin18 (CLDN18) (NM_016369) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Claudin18
Synonyms:	SFTA5; SFTPJ
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC212135 representing NM_016369 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCCACCACCACATGCCAAGTGGTGGCGTTCCTCCTGTCCATCCTGGGGCTGGCCGGCTGCATCGCGG
CCACCGGGATGGACATGTGGAGCACCCAGGACCTGTACGACAACCCCGTCACCTCCGTGTTCCAGTACGA
AGGGCTCTGGAGGAGCTGCGTGAGGCAGAGTTCAGGCTTCACCGAATGCAGGCCCTATTTACCATCCTG
GGACTTCCAGCCATGCTGCAGGCAGTGCAGCCCTGATGATCGTAGGCATCGTCTGGGTGCCATTGGCC
TCCTGGTATCCATCTTTGCCCTGAAATGCATCCGATTGGCAGCATGGAGGACTCTGCCAAGCCAACAT
GACACTGACCTCCGGGATCATGTTTCATTGTCTCAGGCTTTGTGCAATTGCTGGAGTGTCTGTGTTGCC
AACATGCTGGTGACTAACTTCTGGATGTCCACAGCTAACATGTACACCGCATGGGTGGGATGGTGCAGA
CTGTTCCAGACCAGGTACACATTTGGTGC GGCTCTGTTTCGTGGGCTGGGTGCTGGAGGCCTCACACTAAT
TGGGGGTGTGATGATGTGCATCGCTGCCGGGGCTGGCACCAGAAGAAACCAACTACAAAGCCGTTTCT
TATCATGCCTCAGGCCACAGTGTTCCTACAAGCCTGGAGGCTTCAAGGCCAGCACTGGCTTTGGGTCCA
ACACCAAAAACAAGAAGATATACGATGGAGGTGCCCGCACAGAGGACGAGGTACAATCTTATCCTTCCAA
GCACGACTATGTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RC212135 representing NM_016369
 Red=Cloning site Green=Tags(s)

MSTTTTCQVVFLLSILGLAGCIAATGMDMWSTQDLYDNPVTSVVFQYEGLWRSCVRQSSGFTECRPYFTIL
 GLPAMLQAVRALMIVGIVLGAIGLLVSIFALKCIRIGSMEDSAKANMTLTSGIMFIVSGLCAIAGVSVFA
 NMLVTNFWMSTANMYTGMGGMVQTVQTRYTFGAALFVGWVAGGLTLIGGVMMCIACRGLAPEETNYKAVS
 YHASGHSVAYKPGGFKASTGFGSNTKNNKIYDGGARTEDEVQSYPYSKHDYV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_016369

ORF Size: 783 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_016369.4
RefSeq Size:	3359 bp
RefSeq ORF:	786 bp
Locus ID:	51208
UniProt ID:	P56856
Cytogenetics:	3q22.3
Domains:	PMP22_Claudin
Protein Families:	Transmembrane
Protein Pathways:	Cell adhesion molecules (CAMs), Leukocyte transendothelial migration, Tight junction
MW:	27.7 kDa
Gene Summary:	<p>This gene encodes a member of the claudin family. Claudins are integral membrane proteins and components of tight junction strands. Tight junction strands serve as a physical barrier to prevent solutes and water from passing freely through the paracellular space between epithelial or endothelial cell sheets, and also play critical roles in maintaining cell polarity and signal transductions. This gene is upregulated in patients with ulcerative colitis and highly overexpressed in infiltrating ductal adenocarcinomas. PKC/MAPK/AP-1 (protein kinase C/mitogen-activated protein kinase/activator protein-1) dependent pathway regulates the expression of this gene in gastric cells. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jun 2010]</p>

Product images:



Circular map for RC212135