

Product datasheet for RC212177A1V

Human TL1A (TNFSF15) (NM_005118) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TL1A (TNFSF15) (NM_005118) AAV Particle
Tag:	Myc-DDK
Symbol:	TL1A
Synonyms:	TL1; TL1A; TNLG1B; VEGI; VEGI192A
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC212177 representing NM_005118. Blue=ORF Red=Cloning site Green=Tag(s)

GCTCGTTTAGTGAACCGTCAGAATTTTGTAAACGACTCACTATAGGGCGGCCGGAATTCGTCGACTG
 GATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC
 ATGGCCGAGGATCTGGGACTGAGCTTTGGGAAACAGCCAGTGTGGAATGCTGCCAGAGCACGGCAGC
 TGCAGGCCCAAGGCCAGGAGCAGCAGCGCAGCTGGGCTCTCACCTGCTGCCTGGTGTGCTCCCTTC
 CTTGCAGGACTCACCATACCTGCTTGTGAGCCAGCTCCGGGCCAGGGAGAGGCCTGTGTGCAGTTT
 CAGGCTCTAAAAGGACAGGAGTTTGACCTTACATCAGCAAGTTTATGCACCTCTTAGAGCAGACGGA
 GATAAGCCAAGGGCACACCTGACAGTTGTGAGACAACTCCACACAGCACTTTAAAAATCAGTTCCCA
 GCTCTGCACTGGGAACATGAAGTGGCTGGCTTACCAAGAACCAGTGAAGTATACCAACAAATTC
 CTGCTGATCCCGAGTCCGGGAGACTACTTCATTTACTCCAGGTACATTCCGTGGGATGACCTCTGAG
 TGCAGTGAAATCAGACAAGCAGGCCGACCAAAAGCCAGACTCCATCACTGTGGTCATCAACCAAGTA
 ACAGACAGCTACCCTGAGCCAACCCAGCTCCTCATGGGGACCAAGTCTGTATGCGAAGTAGGTAGCAAC
 TGGTTCAGCCCATCTACCTCGGAGCCATGTTCTCCTTGCAAGAAGGGGACAAGCTAATGGTGAACGTC
 AGTGACATCTCTTGGTGGATTACACAAAAGAAGATAAAACCTTCTTTGGAGCCTTCTTACTA
 ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGAT
 TACAAGGATGACGACGATAAGGTTTAAACGGCCGGC

Protein Sequence:	>Peptide sequence encoded by RC212177 Blue=ORF Red=Cloning site Green=Tag(s)
-------------------	---

MAEDLGLSFGETASVEMLPHEGSCRPKARSSARWALTCCLVLLPFLAGLTTYLLVSQLRAQGEACVQF
 QALKGQEFAPSHQVYAPLRADGDKPRAHLTVVRQTPTQHFKNQFPALHWEHELGLAFTKNRMNYTNKF
 LLIPESGDYFIYSQVTFRGMTSECSEIRQAGRPKNPDSITVVITKVTDSEPTQLLMGTGSVCEVGSN
 WFQPIYLGAMFSLQEGDKLMVNVSDISLVDYTKEDKTFFGAFL
 TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species:	Human
----------	-------


[View online »](#)

Serotype:	AAV-2
ACCN:	NM_005118
ORF Size:	753 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_005118.2</u>
RefSeq Size:	2011 bp
RefSeq ORF:	756 bp
Locus ID:	9966
UniProt ID:	<u>O95150</u>
Cytogenetics:	9q32
MW:	28.1 kDa