

## Product datasheet for **RC212193L3V**

### PTP4A3 (NM\_032611) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PTP4A3 (NM_032611) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PTP4A3
Synonyms:	PRL-3; PRL-R; PRL3
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_032611
ORF Size:	519 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212193).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032611.1</a> , <a href="#">NP_116000.1</a>
RefSeq Size:	1396 bp
RefSeq ORF:	522 bp
Locus ID:	11156
UniProt ID:	<a href="#">O75365</a>
Cytogenetics:	8q24.3
Protein Families:	Druggable Genome, Phosphatase
MW:	19.4 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the protein-tyrosine phosphatase family. Protein tyrosine phosphatases are cell signaling molecules that play regulatory roles in a variety of cellular processes. Studies of this class of protein tyrosine phosphatase in mice demonstrates that they are prenylated in vivo, suggesting their association with cell plasma membrane. The encoded protein may enhance cell proliferation, and overexpression of this gene has been implicated in tumor metastasis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013]