

## Product datasheet for **RC212196L3V**

### CRAT (NM\_000755) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CRAT (NM_000755) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CRAT
Synonyms:	CAT; CAT1; NBIA8
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000755
ORF Size:	1878 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC212196).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_000755.2</a>
RefSeq Size:	2804 bp
RefSeq ORF:	1881 bp
Locus ID:	1384
Cytogenetics:	9q34.11
Domains:	Carn_acyltransf
Protein Families:	Druggable Genome
MW:	70.87 kDa



[View online »](#)

**Gene Summary:**

This gene encodes carnitine O-acetyltransferase, a member of the carnitine acyltransferase family and a key metabolic pathway enzyme which plays an important role in energy homeostasis and fat metabolism. This enzyme catalyzes the reversible transfer of acyl groups from an acyl-CoA thioester to carnitine and regulates the ratio of acyl-CoA/CoA. It is found in both the mitochondria and the peroxisome. Alternative splicing results in transcript variants encoding different isoforms that may localize to different subcellular compartments. [provided by RefSeq, Oct 2016]